



BRONCHOSCOPY ACCESSORISE

CDSCO
(GOVT. OF INDIA)
APPROVED
MANUFACTURER



FDA
R.No.3022121578

Tracheal "Y" Stent



Tracheal Stent



Snare




Temponation Balloon



Silicone Stent



Zero Tip Nitinol /
Hexagonal Basket



Forceps



Mucus Specimen Trap
(Capacity: 70ml.)



HOSPI LINE
EQUIPMENTS PVT. LTD.

A total endoscopy solution.....

Website : www.hospilineindia.com.

ABOUT US

Hospi Line Equipments Pvt. Ltd. is a reputed organization, offering total endoscopy solution, manufacturing biggest range of the products under one roof. Our penetration is in around 20 countries, where we are exporting the products for a long period of time.

Our company has been awarded with MD-5, MD-9, ISO 13485:2016, CE certificates with FDA registered, which is evidence to our deep involvement in the manufacturing of the best quality medical products. Company has a big manufacturing unit in a land area 1800 sq.mtrs., in a very clean environment about 35km from Delhi where air pollution is supposed to be nil, apart of it we assemble and pack our products in one lac and 10k level clean room, constructed in 6000 sq. ft. area. We regularly exhibit in various International exhibitions, in which majors are below:

Arab Health-Dubai	Otolaryngology-Dubai
Medica-Germany	Meditex-Bangladesh
Medexpo-Kenya	Medexpo-Saudi Arabia
TIHE-Tashkent, Uzbekistan	KIHE-Almaty, Kazakistan
Medima-Siberia, Russia	Medi-Pharma-Vietam
World Gastro (WSG)-Ethiopia	World Gastro (WSG)-UAE

Apart of it we have displaying/detailing of products at international platform in the competition of big players of this field gives us excellent feeling and satisfaction about our workmanship and we learn more and more.

We have been awarded on national level two times. We also have in house endoscope service centre for flexible endoscope. All brands (video/fibreoptic) which gives complete knowledge of endoscope anatomy to our engineer which make them very accurate about the function of accessories through the scopes. Our complete team is well aware about the functioning and construction of endoscopes, therefore they know well what they are producing.

Apart of UGI, LGI & ERCP products. Our other products range are:

1. Bronchoscopy accessories
2. Urology Disposable /Instruments/Equipments
3. CCD Camera, Light Source (Halogen / Xenon / LED)

In brief the moto of our company is : An economical approach with latest technology...



1000mp/m³ Cleanroom

TBNA / Aspiration Needle


✓ Disposable



(Lock/Unlock Needle)



(Adjustable Needle)



(Pleural Biopsy Needle)

Channel Cleaning / Washing Brush

✓ Disposable
✓ Reusable




Cytology Brush
(Covered / Open)

✓ Disposable



Biopsy Forcep (Hot / Cold)

✓ Disposable
✓ Reusable
✓ Rotatable/
Non-Rotatable



Foreign Body Forceps & Polypectomy Snare (Hot & Cold)

- ✓ Disposable
- ✓ Rotatable
- ✓ Reusable
- ✓ Non-Rotatable



Channel SizeØ	≥2mm
Sheath Length	100cms.

Drive Sheath Kit

- ✓ Disposable




Particulars	Sheath O.D. / Length	For Channel Size
Drive Sheath	1.90mm / 1050mm	≥ 2mm
Biopsy Forcep	1.50mm / 1150mm	≥ 2mm
Cytology Brush	1.40mm / 1150mm	≥ 2mm

For detailed Specification, please contact us.

Needle Knife


✓ Disposable



Heating Coagulation Probe (With & Without Injection Needle)

✓ Disposable

✓ Reusable



Temponation Balloon (Single & Triple Lumen)

✓ Disposable



Mucus Specimen Trap (Capacity :70ml.)

✓ Disposable



EBUS Needle

✓ Disposable



Spray Cannula

✓ Disposable





KOBRA[®]
Pro **HD**

KOBRA Digital CCD Camera
with inbuilt 50 / 100watt LED Light Source

Three in One Unit (inbuilt):


- Monitor
- LED Light Source
- CCD Camera

Features:


- White Balance Single Press
- Feather Touch Control
- Imported Soakable Couple
- Composite Output

Rigid Bronchoscopy Set
(Paediatric & Adult)

Bronchoscope



Bronchoscopy Sheath



Grasping Forcep



Crocodile type



Grasping Forcep



Rat Tooth Type



Biopsy Forcep

Optical Grasping Forcep




Crocodile Type



Foreign Body Forcep

Suction Tubes



#Other parts available on demand.

MICROATH



(Accessories for $\geq 1.2\text{mm}$ Channel Size)



**TBNA / Aspiration
Needle**




**Cytology Brush
(Covered/open)**




**Temponation
Balloon**



**F.B. Dormia
Basket**




**Spray
Catheter**




**F.B. Three / Four
Prong**



**Biopsy / Rubber Jaw
Forcep**




**Channel
Cleaning Brush**




**Cutting
Coagulation Probe**

Self Expendable Metal Stents



Stent Delivery System



Bite Block / Mouth Guard

(With / Without Strap)



(Adult)



(Pediatric)



(Adult with Teether)



(With Airway system)

HOSPI LINE EQUIPMENTS PVT. LTD.



Works: G-439, 440, 441 & 442, UPSIDC Phase-II, M.G. Raod, Dhaulana, Hapur-245 101, Uttar Pradesh (India)



hospiline96@yahoo.com



www.hospilineindia.com



+91-98995 01997 (Customer Care)



+91-93102 21998



#Manufacturer can change design and color without any prior information.