



CDSCO
(GOVT. OF INDIA)
APPROVED
MANUFACTURER



UROLOGY



HOSPILINE EQUIPMENTS PVT. LTD.

A Leader in Manufacturing of Endoscopy Accessories & Equipments

Licence Number: MD-9 : MFG/MD/2018/000082 (MD-5) : MFG/MD/2019/000078

Trade Lic. No. 20B & 21B : UP3720B000051

www.hospilineindia.com



URO DORMIA BASKET



- *Segura/Helical*
- *Zero Tip Nitinol*
- *Fully detachable handle*




DJ STENTS



- *Open/Close End*
- *Single/Multi Coil*
- *Color: Blue/White*

URETHERAL CATHETER



STONE CONE






STONE PUNCH



LITHOCLAST



FORCEP FOR CYSTOSCOPY/URS



PCN Catheter Set





Specification :-

Pixel (H x V)	1280 x 1024
Image Sensor	1/3"
Minimum illumination	0.8Lux F1.2 5600K
Electronic Shutter	1/50 ~ 1/120,000 Sec
White Balance	Single Press
LED Source	50 Watt
Coupler	Imported Soakable "C" Mount
Monitor	15.5" LED
Output	Composite
Unit Weight	Approx. 8kgs
Unit Dimension (L x w x H)	452 x 152 x 342 mm

With Freezing Feature.
SET is Portable Unit.

TURP SET



HOSPILINE

HOSPI LINE EQUIPMENTS PVT. LTD.

Regd. office: - F-15A, 3rd Floor, Jindal Complex, Subhash Chowk, Laxmi Nagar, Delhi-110 092

Works: G-439,440, 441 & 442 UPSIDC, Phase -II, M.G. Road, (Near Coca - Cola), Hapur,- 201 015 (U.P)

Marketing Office: 301-302, 3rd Floor, Angel Mega Mall, Kaushambhi, Ghaziabad - 201010, Uttar Pradesh.

E-Mail : hospiline96@yahoo.com

Tel.: +91-98995 01997